

# SEM-VE

## Application Programming Interface

**v 1.0**

Prepared by  
**Fibics Incorporated**  
556 Booth St. Suite 200  
Ottawa, Ontario  
Canada K1A 0G1  
1-613-860-0861  
info@fibics.com

For  
**Carl ZEISS SMT**  
Version 2010-April-26

**This document describes the SEM-VE Application Programming Interface which allows external programs to communicate with the SEM-VE software.**

**Note:** To make effective use of this interface, the user is required to have a firm knowledge of Microsoft Windows™ programming and the Microsoft Component Object Model (COM). Further information is available at <http://www.microsoft.com/com> and in numerous third party publications.

Due to a policy of continuous improvement, the information in this document is subject to change without notice and should not be construed as a commitment on the part of Fibics Incorporated.

Although every reasonable effort is made to ensure correctness, Fibics Incorporated can accept no responsibility for any errors that may appear in this document.

**Copyright © Fibics Incorporated, Ottawa, Ontario, Canada, 2010  
All Rights Reserved**

The content of this document is the property of Fibics Incorporated and may not be reproduced in any form, or be communicated to a third party without prior written permission of Fibics Incorporated.

## 1. Introduction

This document is intended for third party developers wishing to write applications that communicate with the SEM-VE application to acquire high resolution images. The Application Programming Interface (API) consists of an Type Library, registered as “*NPVE3Z Library*” that provides access to the SEM-VE application.

As a COM Type Library the API can be accessed by various programming languages such as Visual Basic, Visual C++ and various scripting languages (VBScript, JavaScript).

The API exposes a COM IDispatch interface with the following properties and methods. When the COM object is created programmatically, it will automatically start the SEM-VE application if it is not already running. The functions and properties can then be called by the third party application. Upon destruction of the object, the SEM-VE application will be closed if it was started by the third party application.

## 2. Installing the NPVE3Z Library

The NPVE3Z Library is part of the SEM-VE or NPVE application. Before using the API, make sure the SEM-VE application is properly installed and that the system is calibrated. The NPVE3Z Library is automatically registered when the SEM-VE or NPVE applications is started.

## 3. Using the API

### *1. Using the Type Library with Visual Studio 2008 in Visual C#*

Create a new Windows Forms project.

To add a reference to a type library:

From the Project menu, select **References**.  
Select the **COM** tab.  
Select the “**NPVE3Z Library**” type library from the list  
Click **OK**.

In your application, define a variable of type NPVE3Z.IFibicsSEMVE

```
NPVE3Z.IFibicsSEMVE veobj;
```

and at some point during the initialization of the application, assign it as

```
veobj = new NPVE3Z.FibicsSEMVE();
```

You can then use this object to call any of the API functions described in the following section. The following code excerpt shows a simple example of how to create the object and call the AcquireImage function:

```
namespace WindowsFormsApplication1
{
    public partial class Form1 : Form
    {
        NPVE3Z.IFibicsSEMVE veobj;

        public Form1()
        {
            InitializeComponent();
        }

        private void button1_Click(object sender, EventArgs e)
        {
            try
            {
                veobj = new NPVE3Z.FibicsSEMVE();
                button2.Enabled = true;
                button1.Enabled = false;
            }
            catch (Exception ex)
            {
                MessageBox.Show("Fibics SEM-VE Object creation failed: " +
                                ex.Message);
            }
        }

        private void button2_Click(object sender, EventArgs e)
        {
            try
            {
                veobj.AcquireImage(Convert.ToInt32(textBoxWidth.Text),
                                    Convert.ToInt32(textBoxHeight.Text),
                                    Convert.ToSingle(textBoxDwell.Text),
                                    textBoxFileName.Text);

                while (veobj.Busy)
                    System.Threading.Thread.Sleep(200);
                MessageBox.Show("Image complete");
            }
            catch (Exception ex)
            {
                MessageBox.Show("Image acquisition failed: " + ex.Message);
            }
        }
    }
}
```

```
}  
}
```

## 2. Using the Type Library with Delphi (BDS 2006)

Create a new VCL Forms Application.

From the **Component** menu, select **Import Component...**

Select **Import a Type Library**

Select *NPVE3Z Library* from the list

In the **Palette Page** list, type "Fibics", this will cause the **Generate Component Wrappers** to be checked

Select **Create Unit** and press **Finish**.

When this is done, you can select the FibicsSEMVE object from the **Tool Palette** and drop it onto your main form or create an IFibicsSEMVE object programmatically and call the API functions it exposes, as illustrated in the following short example.

```
var  
    SEMVE: IFibicsSEMVE;  
  
procedure TForm1.Button1Click(Sender: TObject);  
begin  
    try  
        SEMVE := CoFibicsSEMVE.Create();  
        Button2.Enabled := True;  
    except  
        on e: Exception do  
            ShowMessage('Fibics SEM-VE Object creation failed: '+e.Message);  
        end;  
    end;  
  
procedure TForm1.Button2Click(Sender: TObject);  
begin  
    try  
        if (SEMVE <> nil) then  
            SEMVE.AcquireImage(StrToInt(WidthEdit.Text),  
                               StrToInt(HeightEdit.Text),  
                               StrToFloat(DwellEdit.Text),  
                               PChar(FileEdit.Text));  
        while (SEMVE.Busy) do  
        begin  
            Application.ProcessMessages();  
            Sleep(200);  
        end;  
        ShowMessage('Image Complete');  
    except  
        on e: Exception do  
            ShowMessage('Image acquisition failed: '+e.Message);  
        end;  
    end;
```

end;

## 4. The API Interface

### Properties

float FOV;

**read** returns the current FOV (field of view) in microns.

**write** sets the current FOV in microns.

bool Busy;

**read** indicates if the scan generator is imaging (TRUE) or Idle (FALSE)

### Methods

```
void AcquireImage([in] long W, H; [in] float Dwell; [in] LPSTR Filename);
```

W and H are the width and height of the image in pixels. Image pixels are square, so this also defines the aspect ratio of the image. For example, 16384 x 16384 results in a square image, but 4000 x 32768 results in a thin vertical image with an aspect ratio of approximately 1:8. The pixel size will be  $FOV/\max(W, H)$ .

Dwell is the dwell time for each pixel, in microseconds.

Filename is the local filename where the image will be stored. The image is stored as a TIFF file with two pages: the first page contains the high resolution image and the second contains a 2048 x 2048 binned low resolution image.

Start acquiring a high resolution image with the given parameters and stores the resulting image in the given location. If `AcquireImage` is called when an image is already being acquired, the current image will end and a new one will be acquired with the new parameters. In general, `Busy` should be polled to determine when the image is complete as the function returns once the imaging starts.

```
void Cancel();
```

Stops scanning immediately. The resulting image is still written to the file but will be completed with blank data.

## Appendix A. NPVE Library Definition in IDL Format

```
[
  uuid(D8FCB34B-2A2F-47E6-8FAE-D2B366076FC7),
  version(1.0),
  helpstring("NPVE3Z Library")
]
library NPVE3Z
{
  importlib("stdole2.tlb");
  [
    uuid(4D3D2A4C-CB90-4A5F-B92A-5A23C24B8E0D),
    version(1.0),
    helpstring("Dispatch interface for FibicsSEMVE Object"),
    dual,
    oleautomation
  ]
  interface IFibicsSEMVE: IDispatch
  {
    [
      propget,      id(0x000000C9)    ]
    HRESULT _stdcall FOV([out, retval] float * Value );
    [
      propput,      id(0x000000C9)    ]
    HRESULT _stdcall FOV([in] float Value );
    [
      id(0x000000CA)    ]
    HRESULT _stdcall AcquireImage([in] long W, [in] long H, [in] float
Dwell, [in] LPSTR Filename );
    [
      propget,      id(0x000000CB)    ]
    HRESULT _stdcall Busy([out, retval] VARIANT_BOOL * Value );
    [
      id(0x000000CC)    ]
    HRESULT _stdcall Cancel( void );
  };

  [
    uuid(405076F3-8B40-45EB-9291-C2692AE983D8),
    version(1.0),
    helpstring("Events interface for FibicsSEMVE Object")
  ]
  dispinterface IFibicsSEMVEEvents
  {
    properties:
    methods:
  };

  [
    uuid(76EAE96-7FCC-4888-928C-F2CC11D84FE9),
    version(1.0),
    helpstring("FibicsSEMVE Object")
  ]
  coclass FibicsSEMVE
  {
    [default] interface IFibicsSEMVE;
    [default, source] dispinterface IFibicsSEMVEEvents;
  };
};
```